

# DAYTIME WEAR RGP LENSES DECISION TREE



Choice of 1 <sup>st</sup> lenses	Spheric		Internal toric		External toric
<b>Géométries</b>	<b>Myopia management</b>				
Aspherical optical zone (toric if internal toric) Adjustable double reservoir periphery <sup>(1)</sup>	<b>PRE AMYOPIC</b>	$\emptyset 10.20 r_0 = K - 0.05$	<b>PRE AMYOPIC TI</b>	$\emptyset 10.20 r_0 = K / r'_0 = K' + 0.15$ or calculated by the laboratory	<b>PRE AMYOPIC TE</b> 1 <sup>st</sup> trial with PRE AMYOPIC
	<b>Myopia – Hypermetropia – Astigmatism</b>				
Spherical optical zone (toric if internal toric) Adjustable double reservoir periphery <sup>(1)</sup>	<b>PRE PERFORM</b>	$\emptyset 10.20 r_0 = K - 0.05$	<b>PRE PERFORM TI</b>	$\emptyset 10.20 r_0 = K / r'_0 = K' + 0.15$ or calculated by the laboratory	<b>PRE PERFORM TE</b> 1 <sup>st</sup> trial with PRE PERFORM
Spherical optical zone (toric if internal toric) Adjustable aspheric periphery <sup>(1)</sup>	PRE AS XL	Calculated by Click & Fit or $\emptyset 10.20 r_0 = K + 0,05$	PRE AS XL TI	Calculated by Click & Fit or calculated by the laboratory or $\emptyset 10.20 r_0 = K / r'_0 = K' + 0.15$	PRE AS XL TE 1 <sup>st</sup> trial with PRE AS XL
Aspheric optic zone Adjustable spheric periphery <sup>(2)</sup>	PRE DS	Calculated by Click & Fit or $\emptyset 9.50 r_0 = Km + 0.10$			PRE DS TE 1 <sup>st</sup> trial with PRE DS
	<b>Presbyopia</b>				
Aspheric optic zone Aspheric clearances	<b>MVB</b>	$\emptyset 9.60 r_0 = Km - 0.10$	<b>MVB TI</b>	$\emptyset 9.60 r_0 = K - 0,30 / r'_0 = k' - 0.15$ or calculated by the laboratory	<b>MVB TE</b> 1 <sup>st</sup> trial with MVB
Spherical optical zone (toric if internal toric) Adjustable spheric periphery <sup>(2)</sup>	<b>EXPERT DS</b>	Calculated by Click & Fit or $\emptyset 9.00 r_0 = Km$	<b>EXPERT TI</b>	Calculated by Click & Fit or calculated by the laboratory or $\emptyset 9.00 r_0 = K - 0,05 / r'_0 = K' + 0.10 + \text{Axe K}$	<b>EXPERT DS TE</b> 1 <sup>st</sup> trial with EXPERT DS

(1) Adjustable eccentricities in steps of 0.10 (2) Adjustable clearances in step of 0,05

# RGP LENSES FITTING EVALUATION

## ■ Visual acuity monocular/binocular

- For the myopia management range:
- > Optimization of the visual acuity should not be considered before 2 weeks of wear
- > In case of insufficient vision, check the pupillary  $\emptyset$

## ■ Slit lamp examination::

- Centring/ mobility
- Integrity of the anterior segment of the cornea
- Fluorescent image after lens recentering

		UNIFOCALS			PRESBYOPIA	
ALIGNED	PRE AMYOPIC/ PERFORM	PRE AS XL	PRE DS	EXPERT DS IF expected on the lower eyelid	MVB	
		PRE AS TI	PRE DS on toric cornea	After recentering		
TIGHT	PRE AMYOPIC/ PERFORM	PRE AS XL	PRE DS	EXPERT DS after recentering	MVB	
FLAT	PRE AMYOPIC/ PERFORM	PRE AS XL	PRE DS	EXPERT DS after recentering	MVB	

Technical support:

Ph. +33 (0)1 45 13 18 45 - Email : [technique@precilens.com](mailto:technique@precilens.com)