

Range	Label reference	Claim	Wearing condition	Renewal	Possible materials
AMYOPIC	AMYOPIC ENFANT	Myopia control	Daytime	3 months	Definitive 74% / GM Advance 49%
AMYOPIC	AMYOPIC TORIC ENFANT	Myopia control	Daytime	3 months	Definitive 74% / GM Advance 49%
AMYOPIC	AMYOPIC JEUNE	Myopia control	Daytime	3 months	Definitive 74% / GM Advance 49%
AMYOPIC	AMYOPIC TORIC JEUNE	Myopia control	Daytime	3 months	Definitive 74% / GM Advance 49%
AURA	AURA	Ametropia	Daytime	12 months	GM Advance 49%
AURA	AURA TS	Ametropia	Daytime	12 months	GM Advance 49%
AURA	AURA TD	Ametropia	Daytime	12 months	GM Advance 49%
C2	C2TT	Ametropia	Daytime	3 months	Definitive 74% / GM Advance 49%
C2	C2NVS	Ametropia	Daytime	12 months	GM Advance 49%
C2	C2NVST	Ametropia	Daytime	12 months	GM Advance 49%
C2	C2T	Ametropia	Daytime	3 months	Definitive 74% / GM Advance 49%
DIVA	DIVA	Ametropia	Daytime	12 months	GM3 58%
DIVA	DIVA TN	Ametropia	Daytime	12 months	GM3 58%
EVA	EVA MULTIFOCAL DISTANCE	Ametropia	Daytime	3 months	Definitive 74% / GM Advance 49%
EVA	EVA MULTIFOCAL NEAR	Ametropia	Daytime	3 months	Definitive 74% / GM Advance 49%
EVA	EVA MULTIFOCAL TORIC DISTANCE	Ametropia	Daytime	3 months	Definitive 74% / GM Advance 49%
EVA	EVA MULTIFOCAL TORIC NEAR	Ametropia	Daytime	3 months	Definitive 74% / GM Advance 49%
EVA	EVA SPHERIC	Ametropia	Daytime	3 months	Definitive 74% / GM Advance 49%
EVA	EVA TORIC	Ametropia	Daytime	3 months	Definitive 74% / GM Advance 49%
LECTRA	LECTRA VL CENTRALE	Ametropia	Daytime	12 months	GM Advance 49%
LECTRA	LECTRA VP CENTRALE	Ametropia	Daytime	12 months	GM Advance 49%
LECTRA	LECTRA TORIQUE VL CENTRALE	Ametropia	Daytime	12 months	GM Advance 49%
LECTRA	LECTRA TORIQUE VP CENTRALE	Ametropia	Daytime	12 months	GM Advance 49%
M44	M44 M1	Ametropia	Daytime	12 months	GM Advance 49% / Contaflex 38%
M44	M44 M2	Ametropia	Daytime	12 months	GM Advance 49% / Contaflex 38%
M44	M44 M3	Ametropia	Daytime	12 months	GM Advance 49% / Contaflex 38%
MEDIGEL	MEDIGEL	Ametropia	Daytime	12 months	GM Advance 49% / Contaflex 38%
MINIGEL	MINIGEL	Ametropia	Daytime	12 months	GM Advance 49% / Contaflex 38%
SB60	SB60	Ametropia	Daytime	12 months	GM3 58%
VERFLEX	VERFLEX SPHERIQUE	Ametropia	Daytime	12 months	GM Advance 49% / Contaflex 38%
VERFLEX	VERFLEX C	Ametropia	Daytime	12 months	GM Advance 49%